

Restoring Focalpoint® 3.0 Window Settings.

From within Focalpoint®, click **File**, then click **Open**.

Under Directories (Windows 3.x) or Folders (Windows 95), proceed according to what is displayed on the screen:

If there is a folder named DATADIR, double-click the folder.

If no DATADIR folder is visible, double-click the driver letter G:

Double-click the **APOWSTN** folder.

Under File Name, double-click **FPOINT.FPI**.

(Wait as Focalpoint® reloads the original default configuration)

From the File menu, click **Save As**. Then click **OK**.

Click **OK** to confirm overwriting the **FPOINT.FPI** file.

The original settings have now been restored to your workstation and will remain in effect when Focalpoint is restarted.